

MULTI-MODALITY

**COMPLETE
NEURODIAGNOSTIC
Solutions
EEG, EMG, NCS, EP
& ULTRASOUND**




Cadwell neurodiagnostic solutions are engineered to work on a single cart or as separate systems connected with CadLink[®] data management.

Cadwell is dedicated to engineering smart and affordable neurodiagnostic technology.

ONE SYSTEM THAT CAN DO IT ALL

Tailor-made



Electromyography	EMG Single Fiber EMG Multi Motor Unit Analysis Turns / Amplitude Analysis	Interference Pattern Analysis Macro EMG Guided Injection Protocol
Nerve Conduction	Motor / Sensory / Mixed / CSI F Wave / H Reflex Repetitive Nerve Stimulation	Blink Reflex MUNE
Evoked Potentials	SEPs Dermatomal SEPs Interleave SEPs VEP/ ERG BAER	MLR LLR P300 VEMP ECochG
Advanced Protocols	Collisions Studies Custom Protocols	Trigger In x 4 Trigger Out x 4
Autonomic Studies	Heart Rate Variability Studies	Sympathetic Skin Response
Electroencephalography	Routine EEG Ambulatory EEG	cEEG for ICU EMU / LTM
Ultrasound Imaging Modes	B, 2B, 4B, BM, M Zoom, B-Steer Compound Trapezoid Color Doppler	Power / Directional Power Spectral Doppler Duplex, Triplex HPRF, THI, CW Auto Image Optimization
Ultrasound Transducers	Linear Phased Array	5.0 - 18.0 MHz
Configurations	T2 Cart - single or dual monitor Point of Care Cart	Portable wheeled case Portable carry case